

INJECTABLE ROOT CANAL SEALER

R-ROOT SP

VETERINARY

KEY FEATURES

- PREMIXED FORMULA FOR EASY APPLICATION
- MOISTURE-ACTIVATED SETTING
- SLIGHT EXPANSION UPON CURING FOR AN EXCELLENT SEAL
- HIGH INITIAL ALKALINITY WITH ANTIBACTERIAL PROPERTIES





INDICATIONS

- Permanent obturation after root canal preparation of vital teeth
- Permanent obturation after root canal preparation of non-vital or infected teeth
- Compatible with Lateral Condensation, Single Cone, and Warm Vertical Compaction techniques

DIRECTIONS FOR USE

- Step 1** Clean and thoroughly dry the root canal.
- Step 2** Attach the dispensing tip to the syringe.
- Step 3** Insert the tip into the canal and dispense the filling material.