


B E T T E R
O R A L
B E T T E R
L I F E

 BETTER DENTAL BETTER LIFE®
奥凯伟迪



CONTENTS

4	About OKVD
5	Infinite Bond
5	Instant Bond Universal
6	Flow Fill
6	R-Root SP
7	Easy Fix Periodontal Adhesive
7	EasyOrtho Orthodontic Adhesive
8	Blue Glue
8	Clear Aligner Attachment Bonding Adhesive
9	Addition Silicone (VPS) Impression Material
9	Putty with Light Body Bundle (Impression Material)
10	Putty (Impression Material)
10	Heavy Body Impression Material
11	Heavy Body Elastomeric Impression Material
11	Light Body Impression Material
12	Transparent Silicone
12	Transparent Silicone via Dispenser Gun
13	Bite Registration
13	Gum reproduction Silicone
14	Silicone Mixer Machine
14	Silicone Dispensing Gun
14	Tray Adhesive
15	Transparent Tray with Elastic Film for Dental Implant
16	Transparent Tray (Orthodontics + Aesthetic Restoration)
16	Rubber Dam Sheet
17	Rubber Dam Toolkit

ABOUT OKVD

Beijing OKVD Biological Technology Co., Ltd was founded in 2011, located in Beijing airport economic core area. OKVD is GMP certified and has passed ISO13485:2016. OKVD focuses on R&D and manufacturing of dental products. With an excellent R&D team including dozens of outstanding experts and professors, OKVD is dedicated to become a global reputable, reliable dental products supplier.

- No 1. Rubber Dam product manufacturer in China
- One of the largest silicone impression material manufacturer in China
- Approved Supplier for No.1 Dental hospital in China (Peking University School and Hospital of Stomatology)
- Trusted by more than 8000 dental clinics and hospitals.
- Provide private label service internationally



Infinite Bond



Main Features :

- Moderate flexibility;
- Organic boron polymerization catalyst, suitable for oral environment bonding;
- Low monomer residue, producing excellent biocompatibility;
- Strong bonding performance with various dental materials, including titanium alloy, PEEK, ceramic, and zirconia .etc.

Clinical Applications :

- Stabilization of loose teeth;
- Bonding for inlays, crowns, and bridges;
- Root canal fillings and posts bonding.

Order Info

Product Name	Size	Specs	Order No
Infinite Bond	PC	1 X Catalyst 0.7g 1 X Polymer (Clear) 3.0g 1 X Polymer (Radiopaque) 5.0g 1 X Teeth Primer 4.0g 2 X Monomer 5.0g / Pc 1 X Dispensing Dish 2 X Dispensing Cups 1 X Syringe 1 X Dental Needles 3 X Needle Sealing Stopper 1 X Measuring Spoon 1 X Brush Handle 10 X Brush Tips (Pink.Brush-Dip) 10 X Brush Tips (Blue.Bulk-Mix)	75001

Instant Bond Universal



Main Features :

- One-step light-curing reduces operation time;
- Excellent bonding performance.

Clinical Applications :

- Bonding for composite resin restorations;
- Core buildup for abutment tooth reconstruction;
- Surface treatment for dental posts and other dental materials.

Order Info

Product Name	Size	Specs	Order No
Instant Bond Universal	PC	5ml/bottle	55001

Flow Fill



Main Features:

- Excellent flowability, ensuring void-free obturation;
- High mechanical strength, tough and wear-resistant;
- Moderate radiopacity.

Clinical Applications :

- Filling of micro-cavities;
- Contouring of complex cavities;
- Base filling for large cavities.

Order Info

Product Name	Size	Specs	Order No
Flow Fill flowable resin	PC	2.2g/syringe, 5 x tips	45001

R-Root SP



Main Features :

- Premixed, hydraulic setting, easy to use;
- Slight expansion upon setting, ensuring superior sealing;
- Strong alkalinity before curing, providing antibacterial effects.

Clinical Applications :

- Permanent filling after vital pulp root canal preparation;
- Permanent filling after necrotic teeth or infected root canal preparation;
- Suitable for lateral condensation, single-cone, and vertical compaction techniques.

Order Info

Product Name	Size	Specs	Order No
R-Root SP	PC	1g/syringe, 6 x tips	65001

| Easy Fix Periodontal Adhesive



Main Features:

- Easy operation, short setting time (about 20 seconds)
- Excellent flexibility. Provide a certain degree of movement, cushioning, and absorption of deformation stress for loose teeth.
- Light cure
- Easy to maintain oral hygiene environment

Order Info

Product Name	Specs	Order No
Easy Fix Periodontal Adhesive	1x2g paste, 6xtips, 1xprimer	035001

| EasyOrtho Orthodontic Adhesive.....



Main Features:

- No primer needed, reduce operational errors and save clinical operation time
- Excellent bonding performance
- Used on all type of bracket (metal , plastic, ceramic etc.)
- Mainly composed of methacrylic acid monomer, glass powder, photopolymerization initiator
- Light cure

Order Info

Product Name	Size	Specs	Order No
EasyOrtho Orthodontic Adhesive	PC	4g paste (primer included)	17004
EasyOrtho Orthodontic Adhesive	PC	2g paste (primer included)	17005

Blue Glue



Main Features:

1. Specialized Orthodontic Bite-Raising Resin

Specifically designed for orthodontic treatment to increase molar crown height and facilitate bite opening.

2. Non-Stick Handling

Exhibits excellent non-stick properties, minimizing instrument adhesion and improving clinical operability.

3. High Sculptability

Provides easy manipulation and precise contouring for accurate occlusal build-up.

4. High Bond Strength

Delivers reliable adhesion to tooth surfaces, ensuring stability during orthodontic treatment.

Clinical Applications

- For bite opening and orthodontic band cementation.

Order Info

Product Name	Size	Specs	Order No
Blue Glue	PC	5g/syringe	017006

Clear Aligner Attachment Bonding Adhesive



Main Features:

1. Resin Specifically Designed for Clear Aligner Attachments

Developed for the fabrication and bonding of attachments used in clear aligner orthodontic treatment.

2. Non-Stick Handling Properties

Demonstrates excellent non-stick behavior, reducing instrument adhesion and improving clinical handling efficiency.

3. Easy Demolding After Curing

Allows clean and effortless removal of attachment templates following polymerization, ensuring precise attachment morphology.

4. Reliable Bond Strength with Primer

When used in combination with a primer, provides dependable adhesion to enamel surfaces.

Clinical Applications

- For enamel bonding of orthodontic attachments used with clear aligners.

Order Info

Product Name	Size	Specs	Order No
Clear Aligner Retainer Attachment Adhesive	SET	4g/syringe+2g/bottle	017007

Addition Silicone (VPS) Impression Material

Putty with Light Body Bundle (Impression Material)



Main Features:

- it is suitable for two-phase primary impression method and two-phase secondary impression making technology;
- It ensures clear and highly accurate impressions.

Order Info

Product Name	Specs	Order No
Putty Set	1x400g Base, 1x400g Catalyst, 1x50ml Light Body, 2xSpoons, 6xMixing tips, 6xNeedle tips	201014

Putty

Impression Material



Order Info

Product Name	Specs	Order No
Putty	1x400g Base, 1x400g Catalyst, 2 Spoons	201015

Main Features:

- With the use of detailed impression material, it is suitable for two-phase primary impression method and two-phase secondary impression making technology, clear and highly accurate;
- The material has good elasticity, and easy to be taken out from the mouth;
- Excellent volume stability, the impression will not be deformed for 1 year, can be repeated to fill the mold;
- Impression can be cleaned, soaked and sterilized under high temperature and high pressure;
- Odorless and tasteless, no mucosal irritation, good biocompatibility.

Heavy Body

Impression Material



Main Features:

- Accurate dose, high material mixing quality, uniform and no bubbles;
- Operation sensitivity is very low, no worry about latex gloves and other pollution;
- Strong hydrophilic, with a high affinity for wet teeth and oral soft tissue;
- High deformation recovery rate, good volume stability, no deformation in 1 year, can repeat the filling mold;
- Can be cleaned, soaked and pasteurized at high temperature.

Order Info

Product Name	Size	Specs	Order No
Heavy Body	380ml/ea	2 X Cartridges 1 Locking Ring	201012

Heavy Body

Elastomeric Impression Material



- Main Features:**
- Easy, efficient handling with precise dosing;
 - Bubble-free mixing and excellent detail reproduction;
 - Strong hydrophilicity for reliable results on wet surfaces;
 - High tear strength, tensile strength, and deformation recovery;
 - Long-term dimensional stability, accurate even after 1 year;
 - Safe with gloves, resistant to contamination;
 - Sterilizable by cleaning, soaking, or high-temperature pasteurization.

Order Info			
Product Name	Size	Specs	Order No
Heavy Body	50ml	2 x Cartridges	201013-1

Light Body

Impression Material



- Main Features:**
- Can be divided into conventional flow, High flow and ultra high flow;
 - Hydrophilic, with wet teeth and oral soft tissue;
 - A high degree of affinity;
 - Excellent fluidity, can penetrate into the fine structure of teeth.

Order Info			
Product Name	Size	Specs	Order no
Light Body Super High Flow (V-light)	50ml/ea	2 X Light body + 2X Mixing Tips +	201007
Light Body High Flow (Light)	50ml/ea	12X Oral Needle Tips	201009
Light Body High Flow (Fast)	50ml/ea		201008

Transparent Silicone.....



Main Features:

- The world’s first machine-mixed transparent silicone;
- Designed for implant repair, transparent and visible, used for transfer rod, abutment transfer and implant location determination;
- With the transparent tray, the use effect is better, easy to judge the implant position, can make the complex implant repair molding procedure become easy;
- Can be used for making various guide plates in dental aesthetic restoration;
- Excellent elasticity, easy oral extraction;
- Gypsum can be easily separated from transparent silicone.

Clinical Applications:

- Implant restoration;
- Applications for dispenser gun mixed transparent silicone.

Order Info

Product Name	Size	Specs	Order No.
Transparent Silicone	380ml	1xCatridge	201011-1

Transparent Silicone via Dispenser Gun



Main Features:

- Non-obscurant, class A addition molding silicone;
- Strong hydrophilic, with a high affinity to wet teeth and oral soft tissue;
- Injection style, easy to operate in the mouth;
- Sterilizable.

Clinical Applications:

- Guides and carriers for aesthetic restoration of dental defects;
- Seal denture;
- X-ray ball fixation during implant depth design;
- Guides for bite reconstruction;
- Guides for indirectly bonded orthodontic brackets.

Order Info

Product Name	Size	Specs	Order No
Transparent Silicone via Dispenser Gun	50ml	2 x Cartridges	201010

Bite Registration



Main Features:

- Non-wax material, no muscle deformation will be caused when taking adhesive record;
- Short solidification time <1 min;
- High strength, Shore hardness up to 90 or more;
- There is no deformation after solidification;
- Good volume stability (no sensitivity to heat);
- Sterilizable.

Clinical Applications:

- Record various conventional bonding;
- Record face bow transfer bonding;
- Segment model replication.

Order Info

Product Name	Size	Specs	Order No
Bite Registration	50ml	2 x Cartridges	201006

Gum Reproduction Silicone



Main Features:

- Good elasticity, used for direct or indirect method and reproduction technology, artificial teeth can accurately and aesthetic reflect the shape and position of the dental tissue around the implant neck or abutment;
- High volume stability and good elasticity;
- Strong tear resistance, can be used repeatedly.

Clinical Applications:

- Reproduction of dental heel tissue around implants in working models;
- Reproductions of heel tissue around abutments during crown/bridge restorations.

Order Info

Product Name	Size	Specs	Order No
Bite Registration	50ml	2 x Cartridges	201006
Gum Reproduction	50ml	2 x Cartridges	201005

Silicone Mixer Machine



Main Features:

- 5:1 ratio injection equipment for machine-mixed silicone;
- Aviation aluminum alloy shell, high structural strength;
- Highly integrated line control;
- High quality mixing and easy cleaning;
- Simple and clear function buttons, one-key start, one-key return.

Specs		Specs	
Size	46*27*16	A-Weighting Sound	<= 72dB(A)
Weight	6.7kg	Pressure Level	
Temperature	0-45	Mixing Time	16 sec
Voltage	AC 100V-240V	Top Thrust	500kg
Frequency	50Hz-60Hz	Mixing Speed	2 speed

Silicone Dispensing Gun



Main usage:

1:1 or 2:1 type rubber tube type silicone manual dispensing equipment.

Order Info

Product Name	Size	Specs	Order No
Silicone Mixer	PC	1 pc/box	017002
Silicone Dispenser Gun	PC	1 pc/box	018001

Tray Adhesive



Main Features:


- Can effectively improve the silicone strength between impression material and tray (metal and plastic tray), thus improving the accuracy of impression;
- Prevent the sliding and falling off of silicone material on the tray;
- Simple coating and convenient operation.

Order Info

Product Name	Size	Specs	Order No
Tray Adhesive	10ml	1 bottle/box	016001

Transparent Tray

with Elastic Film



Main Features:

- Hard plastic design, high strength, small deformation;
- Transparent film on top, visibility is strong, can be penetrated, easy for implant repair transfer studs during implant repair out, ensure the accuracy of impression position;
- Perforated edge design for increased retention of silicone;
- Used for implant restoration impressions, especially for open window impressions.

Use of Open-tray Impression for Implants:

- Place transparent silicone (recommended) or putty evenly on transparent tray;
- Evenly inject light silicone into the oral implant area;
- Place the transparent tray in the entrance and slowly position it. When the implant transfer rod is pushed to the transparent membrane of the tray, use a sharp instrument (probe, etc.) to move the tray. The contact area of the transparent film is punctured. After the transfer rod is pierced, the tray is fixed stably and set for 3 minutes;
- After the silicone in the tray is completely solidified, take it out from the mouth after loosening the transfer rod with special technology;
- Continue to connect the transfer rod to the replacement body and apply silicone coating to the tooth heel. Plaster model can be poured in 1 hour later.

Order Info

Product Name	Size	Specs	Order No
Transparent Tray (Upper Jaw) Extra Large	PC	1/bag	015001
Transparent Tray (Upper Jaw) Large	PC	1/bag	015002
Transparent Tray (Upper Jaw) Medium	PC	1/bag	015003
Transparent Tray (Upper Jaw) Small	PC	1/bag	015004
Transparent Tray (Lower Jaw) Extra Large	PC	1/bag	015005
Transparent Tray (Lower Jaw) Large	PC	1/bag	015006
Transparent Tray (Lower Jaw) Medium	PC	1/bag	015007
Transparent Tray (Lower Jaw) Small	PC	1/bag	015008

Transparent Tray

(Orthodontics + Aesthetic Restoration)



Main Features:

- Fully transparent design, high strength;
- Make silicone more capable of bearing pressure without deformation;
- High light transmittance, does not affect the irradiation of the curing lamp;
- Used for reverse injection stratification technique of flow composite resin in front and back teeth.

Order Info

Product Name	Size	Specs	Order No
Transparent Tray Medium	PC	1/bag	020001

Rubber Dam Sheet



Main Features:

- Three types of materials: ultra-thin, medium and thick;
- Ultra-thin type have good elasticity;
- Medium type with excellent elasticity and tear strength;
- Thick type with super sealing;
- Black blanket is designed for aesthetic restoration.

Order Info

Product name	Specs	Order No
Rubber Dam Sheet 6" x 6"	Pink Medium Thick (peach) 1/bag	104051
	Light Blue Medium Thick (mint) 1/bag	104055
	Purple Thick (Taro) 1/bag	104055
	Green Thick (Apple) 1/bag	104058
	Orange Thick (Orange) 1/bag	104057
	Black Thick (Chocolate) 1/bag	104056
Refill (6" x 6")	Green (thin) 36/box	104004
	Light blue (Medium Thick) 36/box	104008
	Pink (Medium Thick) 36/box	104010
	Black (Medium Thick) 36/box	104012
	Purple (Medium Thick) 36/box	104011
	Green (thin) 52/box	104015
Refill (5" x 5")	Light Blue (Medium Thick) 52/box	104027
	Pink (Medium Thick) 52/box	104028

Rubber Dam Toolkit



Permanent Teeth
Classic Set



Permanent Teeth
Deluxe Set



Deciduous Teeth
Standard Set

Order Info

Product Name	Size	Specs	Order No
Rubber Dam Toolkit classic set (permanent teeth)	set	1 x punch plier + 1 x forceps + 12 x Clamps + 1 x frame + 1 x template + 1 x manual	114001
Rubber Dam Toolkit deluxe set (permanent teeth)	set	1 x punch plier + 1 x forceps + 12 x Blue Clamps + 1 x frame (L) + 1 x frame (M) + 1 x template + 1 x manual + 1 box of Blue Sheet (6" x 6")	114041
Rubber Dam Toolkit Standard set (deciduous teeth)	set	1 x punch plier + 1 x forceps + 11 x Clamps + 1 x frame (S) + 1 x template + 1 x manual	114003
punch plier (5 holes)	PC	1/bag	109001
punch plier (6 holes)	PC	1/bag	109003
Forceps (Straight Head)	PC	1/bag	107001
Forceps (round end)	PC	1/bag	107003
Frame (L)	PC	1/bag	110001
Frame (M)	PC	1/bag	110003
Frame (S)	PC	1/bag	110005

Permanent Teeth Classic Set



Type	Clinical Use	Order No
#9	Butterfly clamp, used for Upper and Lower Jaw front Teeth	108018
#0	Premolars	108022
#1	Upper Jaw Premolars	108023
#2	Lower Jaw Premolars	108025
#2A	Premolars. Can be used for Lower Teeth Gum	108026
#7	Lower Molar	108057
#8A	molar	108032
#12A	Upper Left/Lower Right Molar	108060
#13A	Upper Right/Lower Left Molar	108061
#14A	Molar	108038
#26, #28	Back Teeth	108062, 108063

Deciduous Teeth Standard Set



Type	Use	Order No
#0	Front Teeth (Upper)	108022
#00	Front Teeth (Lower)	108016
#54	Upper Right First Deciduous Teeth	108343
#64	Upper Left First Deciduous Teeth	108345
#55	Upper Right First Deciduous Teeth	108344
#65	Upper Left Second Deciduous Teeth	108346
#74	Lower Left First Deciduous Teeth	108347
#84	Lower Right First Deciduous Teeth	108349
#75	Lower Left Second Deciduous Teeth	108348
#85	Lower Right Second Deciduous Teeth	108350
#U67	Deciduous Teeth	108351

Clamps for Root Canal (Power set)



Root Canal 10 pcs set
9T,1T,1AT 2T,2AT,12A,13A,14T,56T,8AD

Order Info

Product Name	Specs	Order No
Clamps for Root Canal (Power set)	10 pcs	114070

Clamps for Aesthetic Restoration



Aesthetic Restoration 10 pcs set
 4 x #B4 and 4 x #44. 212L+212R, 2. 2A, 206, 207

Order Info		
Product Name	Specs	Order No
Clamps for Aesthetic Restoration	10 pcs	114069

Permanent Teeth refill 8 pcs set



Type		Order No
#212	Butterfly clamp for Front Teeth, Wedge-Shaped	108020
#1, #2A	Premolars	108023, 108026
#7A, #8, #201, #14	Back Teeth	108058, 108049, 108053, 108037
#18	Clamps without Wings	108040
#B4 Matte, Gold, Blue, Black	Multipurpose for Veneer Adhesion, Wedge-Shaped Defect Restoration.	108070
#212L, #212R	Modified #212 Half Butterfly clamp for Front Teeth	108336, 108337

